



Folding guide

DryMax Super Superabsorbent Dressings


In-between toes / fingers



1 Fold 11x20cm dressing in half lengthways.

2 Weave the dressing in between toes or fingers.

Heel / elbow



1 Fold 11x20cm dressing in half widthways.

2 Fold two top corners inwards. Place tape to hold folded corners in place.



3 Open out folded dressing.

4 Place dressing on heel or elbow.

Leg



- 1** Fold 37x56cm dressing from both sides at the malleolus to cover the foot.



- 2** Continue to fold the dressing around the entire leg.



- 3** Make sure the wound area is covered. Use tape to hold the folded boot shape in place.



- 4** Pass a tube gauze as a sock over the dressing.



- 5** The dressing is in place.

Arm



- 1** Fold 37x56 cm dressing around the arm so that the wound area is covered.



- 2** Thread on a tube gauze over the arm.



- 3** The dressing is in place.

Chest



1 Place the dressing over the chest so that the wound area is covered. If necessary, the dressing can be folded and adjusted to the contours of the body.

2 Make one cut on each side of a tube gauze to make holes for the arms. Thread on the tube gauze to hold the dressing in place. The dressing is in place.
***Only cut in the gauze tube, not in the dressings.**

Abdomen



1 Place the dressing over the abdomen so that the wound area is covered. If necessary, the dressing can be folded and adjusted to the contours of the body.

2 Make one cut on each side of a tube gauze to make holes for the arms. Thread on the tube gauze to hold the dressing in place. The dressing is in place.
***Only cut in the gauze tube, not in the dressings.**

Back



1 Place the dressing over the back so that the wound area is covered. If necessary, the dressing can be folded and adjusted to the contours of the body.

2 Make one cut on each side of a tube gauze to make holes for the arms. Thread on the tube gauze to hold the dressing in place. The dressing is in place.
***Only cut in the gauze tube, not in the dressings.**

Groin, thigh and hip



1 Place the dressing over the thigh or hip so that the wound area is covered. If necessary, the dressing can be folded and adjusted to the contours of the body.

2 Make one cut of a tube gauze to make hole for the leg. Thread on the tube gauze to hold the dressing in place. The dressing is in place.

**Only cut in the gauze tube, not in the dressings.*

Ordering information

Item	Product size cm	REF	Pcs/box	PIP code	NHS cat number	Absorption capacity*
DryMax Super	11x10	F60073/10	10	421-3369	EJA293	130 ml
DryMax Super	11x20	F60075/10	10	421-3377	EJA288	330 ml
DryMax Super	20x20	F60076/10	10	421-3385	EJA285	600 ml
DryMax Super	20x30	F60077/10	10	421-3393	EJA294	960 ml
DryMax Super	37x56	F60079/10	10	421-3401	EJA283	3000 ml

*According to international standard EN 13726

WE'RE HERE TO HELP

Contact our UK partner to learn more about the DryMax superabsorbent range.

Email: info@cdmedical.co.uk

Phone: 01942 816184

FREE SAMPLES!

Experience the difference and compare the outcome with our superabsorbent dressing DryMax Super.

