

Folding guide

DryMax Super a superabsorbent dressing for exuding wounds



In-between toes / fingers




1 Fold 5x5 cm dressing in half.



2 Place the dressing in between toes or fingers.

Heel / elbow



1 Fold 11x20 cm dressing in half widthways.



2 Fold two top corners inwards. Place tape to hold folded corners in place.



3 Open out folded dressing.

4 Place dressing on heel or elbow.

Leg



- 1** Fold 37x56 cm dressing from both sides at the malleolus to cover the foot.



- 2** Continue to fold the dressing around the entire leg.



- 3** Make sure the wound area is covered. Use tape to hold the folded boot shape in place.



- 4** Pass a tube gauze as a sock over the dressing.



- 5** The dressing is in place.

Arm



- 1** Fold 37x56 cm dressing around the arm so that the wound area is covered.



- 2** Thread on a tube gauze over the arm.



- 3** The dressing is in place.

Chest



- 1** Place the 37x56 cm dressing over the chest so that the wound area is covered. If necessary, the dressing can be folded and adjusted to the contours of the body.

- 2** Make one cut on each side of a tube gauze to make holes for the arms. Thread on the tube gauze to hold the dressing in place. The dressing is in place.
***Only cut in the gauze tube, not in the dressings.**

Abdomen



- 1** Place the 37x56 cm dressing over the abdomen so that the wound area is covered. If necessary, the dressing can be folded and adjusted to the contours of the body.

- 2** Make one cut on each side of a tube gauze to make holes for the arms. Thread on the tube gauze to hold the dressing in place. The dressing is in place.
***Only cut in the gauze tube, not in the dressings.**

Back



- 1** Place the 37x56 cm dressing over the back so that the wound area is covered. If necessary, the dressing can be folded and adjusted to the contours of the body.

- 2** Make one cut on each side of a tube gauze to make holes for the arms. Thread on the tube gauze to hold the dressing in place. The dressing is in place.
***Only cut in the gauze tube, not in the dressings.**

Groin, thigh and hip



- 1** Place the 37x56 cm dressing over the thigh or hip so that the wound area is covered. If necessary, the dressing can be folded and adjusted to the contours of the body.



- 2** Make one cut of a tube gauze to make hole for the leg. Thread on the tube gauze to hold the dressing in place. The dressing is in place.

***Only cut in the gauze tube, not in the dressings.**

Ordering information

Item	Product size cm	REF	Pcs/box	PIP code	NHS cat number	Absorption capacity*
DryMax Super	5x5	F60070/10	10	433-8075	ELZ85079	20 ml
DryMax Super	11x10	F60073/10	10	421-3369	EJA293	130 ml
DryMax Super	11x20	F60075/10	10	421-3377	EJA288	330 ml
DryMax Super	20x20	F60076/10	10	421-3385	EJA285	600 ml
DryMax Super	20x30	F60077/10	10	421-3393	EJA294	960 ml
DryMax Super	37x56	F60079/10	10	421-3401	EJA283	3000 ml

*According to international standard EN 13726

WE'RE HERE TO HELP

Contact our UK partner to learn more about the DryMax superabsorbent range.

Email: info@cdmedical.co.uk

Phone: 01942 816184

FREE SAMPLES!

Experience the difference and compare the outcome with our superabsorbent dressing DryMax Super.

SCAN ME!

