Absorbest

DryMax Super

DryMax Super is a sterile, conformable, superabsorbent dressing for exuding wounds. DryMax Super absorbs and retains exudate in the dressing. The exudate is absorbed vertically by the dressing and retained in the core forming a gel, thus minimizing the risk of maceration and maintaining a moist and balanced wound environment.

The seamless cover simplifies positioning over the wound. The design allows the dressing to be applied with either side facing the wound. This feature also makes it possible to stack dressings on each other if extra absorption capacity is required.

The dressing remains dry and comfortable for the patient. Simultaneously it contributes to a moist wound environment. The hydrophilic contact layer transports fluid effectively away from the wound and distributes it evenly into the dressing's core constructed by our DryMax Technology.

The ultra-thin and highly conformable features allows the dressing to be used on all body parts. DryMax Super is designed to work well together with compression therapy.

Application Guide



The dressing must be placed close to the wound, or placed in close combination with a suitable dressing.



dressing with suitable bandage, film or tape.



DryMax Super can be used in conjunction with compression therapy.

Product size cm	Pcs/box	Absorption- capacity*	REF	PIP Code	NHS Cat Number
11x10	10	130 ml	F60073/10	421-3369	EJA293
11x20	10	330 ml	F60075/10	421-3377	EJA288
20x20	10	600 ml	F60076/10	421-3385	EJA285
20x30	10	960 ml	F60077/10	421-3393	EJA294
37x56	10	3000 ml	F60079/10	421-3401	EJA283

*According to international standard EN 13726



For further information, or to request a product sample scan here





- Cost effective
- User friendly
- Increased patient comfort
- Improved wound healing*

*Compared to conventional wound dressings



DryMax[®] technology

Proud to be the UK distributor of DryMax CDMedical 01942 816184 info@cdmedical.co.uk www.cdmedical.co.uk



Manufactured by Absorbest AB Sweden www.absorbest.co.uk